



IHE Integration Statement		Date: Apr. 11 <sup>th</sup> , 2008
Vendor	Product Name	Version
ESAOTE S.p.A.	Artoscan-C / C-scan / E-scan XQ / E-scan Opera / G-scan / S-scan	9.7A
These products are intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Acquisition Modality	Broad Worklist Query
		PPS Exception Management
Patient Information Reconciliation	Acquisition Modality	none
Portable Data for Imaging	Portable Media Creator	Web Content
	Image Display	none
Internet address for vendor's IHE information: <a href="http://www.esaote.com/dicom.htm">http://www.esaote.com/dicom.htm</a>		
Vendor's e-mail address for information about IHE: <a href="mailto:ihe@esaote.com">ihe@esaote.com</a>		
Links to Standards Conformance Statements for the Implementation		
DICOM <sup>®</sup>	<a href="http://www.esaote.com/dicom.htm">http://www.esaote.com/dicom.htm</a>	
Links to general information on IHE		
In North America: <a href="http://www.ihe.net">www.ihe.net</a>		In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>

ESAOTE S.p.A.  
via Siffredi 58  
16153 Genova (Italy)

© Copyright ESAOTE S.p.A. 2008. All rights reserved.

<sup>®</sup> DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.